




HOW TO USE: Print this template on letter paper, landscape. DO NOT SCALE. Print at 100% or disable page scaling. Cut along the dashed lines and remove inner shape. Insert your cat's nose and mouth into the front and side cut-outs to see which mask will fit best.

AeroKat 
Feline Aerosol Chamber



The exclusive **Flow-Vu**[®] Inhalation Indicator moves toward the cat as they inhale and will only move if there is a good seal.



Mask Sizer

Each **AeroKat**[®] Chamber is shipped with a small and medium mask. The mask selected should cover the cat's nose and mouth but not the eyes.

Helping animals around the world breathe easier



Trudell Medical International[®]

www.aerokat.com